

# BVI

## BVI® Safety Knives

World Class Leader in High Quality Surgical Blades

- Single handed activation
- Consistent blade quality
- Safety options include Crescent, Sideport and Slit



# Complete Portfolio




## Phaco Slit/Keratome Knives



These knives are designed to make a precise incision size to accommodate the phaco tip.

**Safety Slit** knives are designed for penetration and accurate incision width.



- Available in a wide range of sizes (1.4 – 3.2 mm)
- Single bevel and double bevel blade configurations
- Reference line guides a multi-planar incision



Single Bevel Incision



Double Bevel Incision




| Knife Sizes     | Reference Line From Tip |
|-----------------|-------------------------|
| 1.8 mm and 2 mm | 1.75 mm                 |
| > 2 mm          | 2 mm                    |



## Safety

- True no-look, single-handed activation
- Integrated retractable shield prevents sharps injuries during handling and use
- Activated with a slider that provides tactile, visual and audible feedback
- Safety shield does not interfere with grip, offering better control during use



## Paracentesis/Sideport Knives

These knives are designed for accurate width stab incisions.

The angled **Safety Sideport** knives create a consistent incision size.

**Optimum® Safety Sideport** knives enable the surgeon to control the incision width by varying the depth of the incision.



## Crescent Knives

These knives are designed to create optimal scleral tunnel architecture. They can also be used for trabeculectomy and lamellar dissections.

Available in angled, bevel up configuration



## Consistency

- Advanced proprietary manufacturing process, delivers consistency in sharpness
- Safety shield protects blade against inadvertent damage ensuring a consistent incision

## Sharpness

- Enables surgeon to maneuver the blade and optimize incision architecture
- Designed to deliver precision

# BVI Safety Knives



## Safety Slit Knives

### Single Bevel

Qty: 10/Bx

| Description   | SKU #  |
|---|--------|
| 1.4 mm Slit, Angled, Single bevel                   | 378214 |
| 1.5 mm Slit, Angled, Single bevel                   | 378215 |
| 1.7 mm Slit, Angled, Single bevel                   | 378217 |
| 1.8 mm Slit, Angled, Single bevel                   | 378218 |
| 2 mm Slit, Angled, Single bevel                     | 378220 |
| 2.2 mm Slit, Angled, Single bevel                   | 378222 |
| 2.4 mm Slit, Angled, Single bevel                   | 378224 |
| 2.5 mm Slit, Angled, Single bevel                   | 370145 |
| 2.6 mm Slit, Angled, Single bevel                   | 370146 |
| 2.75 mm Slit, Angled, Single bevel                  | 378227 |
| 2.8 mm Slit, Angled, Single bevel                   | 378228 |
| 3 mm Slit, Angled, Single bevel                     | 378230 |
| 3.2 mm Slit, Angled, Single bevel                   | 378232 |
| 2.5 mm Kojo Slit, 45 Degrees, Arc height is 0.55 mm | 373025 |
| 3.2 mm Kojo Slit, 45 Degrees, Arc height is 0.71 mm | 373032 |

### Double Bevel

Qty: 10/Bx

| Description                        | SKU #  |
|------------------------------------|--------|
| 2 mm Slit, Angled, Double bevel    | 378820 |
| 2.2 mm Slit, Angled, Double bevel  | 378822 |
| 2.4 mm Slit, Angled, Double bevel  | 378824 |
| 2.5 mm Slit, Angled, Double bevel  | 378825 |
| 2.65 mm Slit, Angled, Double bevel | 378826 |
| 2.75 mm Slit, Angled, Double bevel | 378827 |
| 2.8 mm Slit, Angled, Double bevel  | 378828 |
| 2.85 mm Slit, Angled, Double bevel | 378829 |
| 3 mm Slit, Angled, Double bevel    | 378830 |
| 3.2 mm Slit, Angled, Double bevel  | 378832 |



## Optimum Safety Sideport Knives

Qty: 10/Bx

| Description                    | SKU #  |
|--------------------------------|--------|
| 15 Degree Sideport, Straight   | 378235 |
| 22.5 Degree Sideport, Straight | 378236 |
| 30 Degree Sideport, Straight   | 378237 |
| 45 Degree Sideport, Straight   | 378238 |



## Safety Sideport Knives

Qty: 10/Bx

| Description                             | SKU #  |
|---|--------|
| 1 mm Sideport, Angled 45 degrees        | 378210 |
| 1.2 mm Sideport, Angled 45 degrees      | 378212 |
| 0.6 mm 23G Sideport                     | 378233 |
| 1.15 mm 20G Sideport, Angled 45 degrees | 378231 |



## Safety Crescent Knives

Qty: 10/Bx

| Description                           | SKU #  |
|---------------------------------------|--------|
| 2.5 mm Crescent, 55 Degrees, Bevel up | 378234 |